



Current FDA approach for preclinical vector biodistribution studies



**Steven R. Bauer
CBER, FDA**

Background



- **Origin of concern**
 - **Preclinical data**
 - unexpected persistence in gonads
 - inadequacy of preclinical data
 - **Vector integration**
 - literature examples of plasmid DNA integration

(vector least expected to be of concern)
- **PCR signal in gonads not necessarily in germ cells nor integrated**

Plasmid integration-1



- **B lymphocytes**
 - **In vivo role of B lymphocytes in somatic transgene immunization. Xiong et al 1997. PNAS 94: 6352**
 - **100μg DNA intrasplenic**
 - **persistence 3 months**
 - **integration in B cells (digest:self-ligate:PCR)**

Plasmid Integration-2



- **Lung and kidney**
 - **Safety study and characterization of E1A-liposome complex gene delivery protocol in an ovarian cancer model. Xing et al. 1998. Gene Therapy 5 : 1538**
 - **plasmid:liposome intraperitoneal**
 - **plasmid DNA in lung and kidney 1 1/2 year post administration: not in ovaries**
 - » **integrated (DpnI digestion)**

Assay Purpose and Selection



- **Distribution studies are designed to address two issues:**
 - i) potential for germline alteration
 - ii) potential for toxicity in other tissues and organs
(as for other biological products)
- **Detection method should be suitable for clinical samples as well as in animal samples**

In vivo Treatment



- **male and/or female animals**
 - appropriate to clinical use, potential future clinical use
 - age
 - developmental state
 - group size
- **vector administered by clinical route**
- **dose selection**
 - vehicle control
 - clinically relevant dose
 - exceed clinically relevant dose

Sacrifice Times



- **early - peak vector levels**
- **multiple time points to assess duration of vector persistence**

Tissue Harvest



- **peripheral blood, injection site, gonads (other organs as needed for more general toxicity studies)**
- **harvest using procedures to minimize cross contamination**

Assay Parameters



- **assay significant sample of genomic DNA from each tissue**
- **use spike controls**
- **perform replicate tests**
- **adequate PCR sensitivity**

PCR Sensitivity and Sample Size Influences Risk Assessment



- **Sensitivity**
 - capable of 1 copy/ μ g genomic DNA
- **Sample size**
 - sufficient replicates
- **Statistical Considerations**

PCR product detection methods



- **Southern**
 - shows size and specificity
- **Slot or dot blot**
 - probe indicates specificity
- **Real time**
 - newest technology, kinetic analysis needed
- **All capable of 1 copy/ μ g**

Primer Selection



- **Detection of therapeutic gene not necessary to address germline insertion**
- **Choose unique amplicon**

Minimal PCR Recommendations




- **3 samples/tissue, 1 μg genomic DNA each
(or sufficient replicates to total sampling of 3 μg)**
- **2 samples run unspiked, 1 sample/tissue run
with spiked control**
- **sensitivity <100 copies of vector/1 μg genomic
DNA**
- **assay should be adequate for study of clinical
samples as well as in animal tissues**

Technical Considerations



- **Spiking in gonadal tissue**
 - **sensitivity control**
 - **signal extinction**
 - **contaminants, competition**
 - **purification control**
 - **check for loss of vector during DNA purification**

Signal Extinction



PCR, gel electrophoresis

Impact on Informed Consent



- **statement in informed consent regarding current results, unknown risk of vector dissemination and transmission to germ cells**
- **temporary use of contraception recommended**
- **autopsy requested in treated patients**

Impact of Results on Clinical Development



- **gonadal signal**
 - **if assay is adequate**
 - **not detected at all times**
 - **transiently positive**
 - **persistently positive**

Gonadal Signal Undetectable at all Time Points



- **clinical study may proceed, no restrictions on patient population**

Gonadal Signal Transiently Positive



- **re-evaluate clinical study as to patient population, severity of illness; may proceed if benefit justifies risk**
- **semen analysis in treated males requested in follow-up where applicable/appropriate**

Gonadal Signal Persistently Positive



- **restrict patient population to sterile individuals**
- **semen analysis in treated males requested in follow-up where applicable/appropriate**
- **analyze source of signal**

Persistent Signal Source



- **Contamination**
- **Signal not from contamination**
 - **determine duration of positive signal**
 - **determine cellular source**

Acknowledgements



- **FDA appreciates efforts of the RAC and NIH to support public discussion of these issues and to encourage research to answer these questions**

Primer selection



- **Detection of therapeutic gene not necessary to address germline insertion**
- **Choose unique amplicon**
 - **amplify from gene to vector backbone or vector only**
 - **avoid:**
 - **homology with endogenous gene**
 - **multigene families**

Determination of Signal Source




- **Contamination**
 - **during harvest**
 - re-evaluate procedures for tissue harvest
 - sterile instruments, order of organs procured
 - **from prior PCR**
 - utilize PCR “sterilization” techniques, UDG
 - evaluate with other primer sets

Signal Source




- **Signal not from contamination**
 - **determine duration of positive signal**
 - **determine cellular source**
 - **analyze semen, sperm**
 - **in situ PCR**
 - **other novel techniques and approaches**
 - **platform studies**

Statistical Considerations



Real Time PCR



Vector/DNA Recovery

